


Derungs 
MEDICAL LIGHTING

LIGHTS FOR VETERINARY



DERUNGS® LIGHTING FOR VETERINARY

Derungs® Licht AG

We are a global leading provider of lighting systems and accessory components for medical technology. Our guiding principle is combining innovative lighting technology with functional design.

Derungs® Lighting Designs for the Healthcare Industry and Veterinary

Our future builds on know-how gained from more than 70 years of experience. Based on that experience, motivated employees are providing well developed solutions for appropriate lighting in the healthcare industry. Since 1999 the company has been part of the Waldmann Group. A company with the motto «light is life».

Derungs® Advancement through Innovation

In addition to our standard solutions, we also offer total solutions tailored to the wishes of our customers. Consulting and support are as much a part of our service as comprehensive light planning and individual lighting design. We adapt our product line to the needs of the market through innovative new designs.

Derungs® Standards

- High Light Intensity
- Top Quality
- Modern Design
- High Safety Standards
- Most Current Medical Standards
- Functionality / Service Comfort
- LED Economical / Long Life
- 10 Year Replacement Part Guarantee



LIGHTS FOR VETERINARY

Today, veterinary medicine uses the whole armory of diagnostic methods and therapeutic means. All work done in the examination room calls for an incredibly sharp and concentrated eye, making lighting presently an essential medical assistance in this context. To ensure the alertness and full concentration of any veterinarian, a reliable, shadow-less lighting - without little heat generation - is required. This light must be easy to handle with


one hand. The change and renovation of veterinary clinics give this possibility to better optimize the lighting for the vet, the owner, and even the animal. Many professional lighting solutions exist which create a comfortable and calming environment, to create visual effects in the room, emphasizing the simplicity of modern design and interior arrangements.

Derungs® Product Groups

D^{med}®
Treatment and Examination Lights

D^{care}®
Care and Reading Lights


D^{lite}®
Room Lighting




D^{med}® TREATMENT AND EXAMINATION LIGHTS


LIGHTS FOR VETERINARY

D^{med}® saturn







Duo 65/65 C




65 C



65 W



65 F



Battery 65 F

- About 55 000 lux / 1 m per luminaire head
- Vision not bothered thanks to homogeneous, shadow-free lit area
- Intense light (size of lit area at 1 m distance: ø 16 cm)
- Solid and easy to handle design leaving the patient feeling worry-free
- Sensational energy savings and lighting efficiency
- Simple, effortless positioning of the light with one hand
- Fixed and stable light cone position thanks to its joints system with weight balance

D^{med}® isis





Duo 50/50 C



50 C




50 W




50 F


- About 54 000 lux / 0.8 m per luminaire head
- Light intensity setting with a 2-stages rocker switch (60/100 %)
- An intense lit area, perfect in color rendering, homogeneous, neutral white
- Space-saving, compact and easy-to-use design
- Sensational energy savings and lighting efficiency (small surgical light)
- Fixed and stable light cone position thanks to its joints system with weight balance

D^{med}® halux LED






20 C L1



20 P LX



20 P SX

- 50 000 lux / 0.5 m
- 5 light dimming levels
- Adjustable color temperature (three color shades allow personalized viewing of contrasts)
- Simple control on the head
- No heat (infrared) radiation
- High energy efficiency and luminous efficiency
- Optimized cooling of the lamp head
- LED technology reduces maintenance costs to a minimum while providing high energy efficiency
- Arm system with friction joint for safe and sure positioning of the light

D^{med}® TREATMENT AND EXAMINATION LIGHTS

LIGHTS FOR VETERINARY



50 P SX



50 P FX


- Positioning of the light is risk-free. Even after long usage, the luminaire head is merely warm
- Noticeable energy and cost savings thanks to a reduced electric consumption and a long service life of the light source
- Attractive and sturdy design guaranteeing a long service life
- Extremely easy up-keep. Perfect and simple cleaning thanks to the closed design of the luminaire (closed wiring passage)
- About 48 000 lux / 0.5 m



Two 35/35 P FX



Light head



- Possibility to switch between spot or flood
- 3-stages rocker switch (spot-0-flood)
- Small lights with extremely high light beam
- Switching at will: Between examination light (10°) or reading light (36°)
- About 30 000 lux / 0.5 m (10° Spot)
- About 5 000 lux / 0.5 m (36° Flood)

Production and international sales

Derungs Licht AG
Hofmattstrasse 12
9200 Gossau
Switzerland
Telephone +41 71 388 11 66
Fax +41 71 388 11 77
www.derungslicht.com
mailbox@derungslicht.com

Distribution sales Switzerland

Waldmann Lichttechnik GmbH
Benkenstrasse 57
5024 Küttigen
Switzerland
Telephone +41 62 839 12 12
Fax +41 62 839 12 99
www.waldmann.com
info-ch@waldmann.com